

## Specialized system for organic solvents



Alpha 3-4 LSCbasic



# Product highlights

## Special system Alpha 3-4 LSCbasic



- Coated real-glass chamber**  
borosilicate glass for higher safety
- Flow tube**  
for optimum vapour flow
- Deep-seated ice condenser**  
to improve performance


### Higher safety for the freeze drying of organic, solvent-containing samples

- Unique system set-up especially for organic solvents
- Flexible use for a wide range of different applications with organic solvents thanks to two fixed packages
- Ice condenser temperature of  $-105^{\circ}\text{C}$  with natural refrigerants
- Particularly easy operation with the LSCbasic control system
- Capacitive vacuum sensor
- Controlled vacuum and solvent-resistant vacuum pump
- EPDM seal suitable for solvents
- Integrated valve for inertisation



### Fixed configurations

For detailed information please visit our website:  
[martinchrist.de/alpha-3-4](http://martinchrist.de/alpha-3-4)





**Martin Christ**  
**Gefriertrocknungsanlagen GmbH**

An der Unteren Söse 50  
37520 Osterode am Harz, Germany

Tel. +49 (0) 55 22 50 07-0

Fax +49 (0) 55 22 50 07-12

[info@martinchrist.de](mailto:info@martinchrist.de)

[www.martinchrist.de](http://www.martinchrist.de)